



Magic™ Anti-68kDa Neurofilament polyclonal antibody (DPABH-00070)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Preparation of bovine spinal cord intermediate filaments.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Mouse, Rat, Chicken, Cow, Human
Purification	IgY fraction
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr, WB, ICC/IF
Format	Liquid
Size	100 µl
Buffer	Constituents: Whole serum
Preservative	10mM Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Neurofilaments usually contain three intermediate filament proteins: L, M, and H which are involved in the maintenance of neuronal caliber.
--------------	---

Keywords

NEFL; neurofilament, light polypeptide; NFL; NF-L; NF68; CMT1F; CMT2E; PPP1R110; neurofilament light polypeptide; neurofilament subunit NF-L; neurofilament triplet L protein; neurofilament protein, light chain; neurofilament, light polypeptide 68kDa; ligh
