



# Anti-antithrombin III polyclonal antibody (DPABH-00080)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Bovine Antithrombin III purified from bovine plasma.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Cow
Purification	Protein G purified
Conjugate	Unconjugated
Applications	WB, ELISA, RIA
Format	Liquid
Size	150 µg
Buffer	Constituents: 50% Glycerol, PBS, pH 7.5
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

Introduction	Most important serine protease inhibitor in plasma that regulates the blood coagulation cascade. AT-III inhibits thrombin as well as factors IXa, Xa and XIa. Its inhibitory activity is greatly enhanced in the presence of heparin.
--------------	---

**Keywords**

SERPINC1; serpin peptidase inhibitor, clade C (antithrombin), member 1; AT3; antithrombin-III; ATIII; serpin C1; antithrombin III; serine (or cysteine) proteinase inhibitor, clade C (antithrombin), member 1;

---