



Rabbit Anti-Human DNAJB1 Polyclonal antibody (DPABH-00321)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DNAJB1 fusion protein, sequence: MGKDYYQTLGLARGASDEEIKRAYRRQALRYHPDKNKEPGAEEKFKEIAEAYDVLSDPRK REIFDRYGEEGLKGSGPSGGSGGGANGTSFSYTFHGDPHAMFAEFFGGRNPFDTFFGQRN GEEGMDIDDPFSGFPMGMGGFTNVNFGRSRSAQEPARKKQDPPVTHDLRVSLLEEIYSGCT KKMKISHKRLNPDGKSIRNEDKILTIEVKKGWKEGKITFPKEGDQTSNNIPADIVFLK DKPHNIFKRDGSDVIYPARISLREALCGCTVNVPTLDGRTIPVVFKDVIIRPGMRRKVPGE GLPLPKTPEKRGDLIIIEFVIFPERIPQTSRTVLEQVLPI (1-340 aa encoded by BC019827)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IP, IHC, IF, ELISA
Positive Control	mouse spleen tissue, COLO 320 cells, HEK-293 cells, HeLa cells, human liver tissue, mouse lung tissue
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Preservative	0.02% Sodium Azide

Storage

Store at -20°C. Aliquoting is unnecessary for -20°C storage.

BACKGROUND

Keywords

DNAJB1; DnaJ (Hsp40) homolog, subfamily B, member 1; Hdj1; Sis1; HSPF1; Hsp40; RSPH16B; dnaJ homolog subfamily B member 1; human DnaJ protein 1; dnaJ protein homolog 1; radial spoke 16 homolog B; heat shock 40 kDa protein 1;

GENE INFORMATION

Entrez Gene ID

[3337](#)

UniProt ID

[P25685](#)
