



Anti-FRMD6 (C-terminal) polyclonal antibody (DPABH-00408)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FRMD6 (FERM domain containing 6) is a protein-coding gene. Diseases associated with FRMD6 include alzheimers disease , and squamous cell carcinoma. An important paralog of this gene is FRMD1.
Immunogen	Synthetic peptide corresponding to Human FRMD6 aa 600 to the C-terminus conjugated to Keyhole Limpet Haemocyanin (KLH).Database link: Q96NE9
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituent: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	FRMD6 FERM domain containing 7 [Homo sapiens]
Official Symbol	FRMD6
Synonyms	FRMD6; FERM domain containing 6; EX1; Willin; C14orf31; c14_5320; FERM domain-containing protein 6; expanded homolog; 4.1 ezrin radixin moesin (FERM)-containing protein;
Entrez Gene ID	122786
Protein Refseq	NP_001035946.1
UniProt ID	Q96NE9
Pathway	Hippo signaling pathway.