



Anti-ZDHHC2 (internal region) polyclonal antibody (DPABH-01439)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Palmitoyltransferase specific for GAP43 and DLG4/PSD95.
Immunogen	Synthetic peptide derived from an internal region of Human ZDHHC2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	African green monkey
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituents: 49% PBS, 0.88% Sodium chloride, 50% Glycerol. Note: (without Mg ²⁺ , Ca ²⁺)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Gene Name [ZDHHC2 zinc finger, DHHC-type containing 3 \[Homo sapiens \]](#)

Official Symbol	ZDHHC2
Synonyms	ZDHHC2; zinc finger, DHHC-type containing 2; DHHC2; ZNF372; palmitoyltransferase ZDHHC2; rec; ream; DHHC-2; zinc finger protein 372; reduced expression in cancer protein; zinc finger, DHHC domain containing 2; zinc finger DHHC domain-containing protein 2; reduced expression associated with metastasis protein;
Entrez Gene ID	51201
Protein Refseq	NP_057437.1
UniProt ID	Q9UIJ5
Function	palmitoyltransferase activity; protein-cysteine S-palmitoyltransferase activity; zinc ion binding;