



Anti-IVL (aa 350-450) polyclonal antibody (DPABH-02011)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Part of the insoluble cornified cell envelope (CE) of stratified squamous epithelia.
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 350 - 450 of Human Involucrin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituent: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	IVL involucrin [Homo sapiens]
Official Symbol	IVL
Synonyms	IVL; involucrin;
Entrez Gene ID	3713
Protein Refseq	NP_005538.2
UniProt ID	P07476
Pathway	Corticotropin-releasing hormone;
Function	protein binding, bridging; structural molecule activity;