



Rabbit Anti-Human VGLL4 (a.a. 201-290) Polyclonal antibody (DPABH-02276)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	VGLL4
Immunogen	KLH conjugated synthetic peptide derived from human VGLL4: 201-290/290
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human, Chicken, Dog, Pig, Cow, Horse, Rabbit, Sheep
Purification	affinity purified by Protein A
Conjugate	Unconjugated
Applications	WB WB: 1:500-2000
Format	Liquid
Concentration	1 mg/ml
Size	50 µl, 100 µl, 200 µl
Buffer	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Preservative	0.02% Proclin300
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

BACKGROUND

Introduction Predicted to enable transcription coactivator binding activity. Involved in negative regulation of Wnt signaling pathway; negative regulation of cell growth; and negative regulation of hippo signaling. Predicted to be located in nucleus.

Keywords Transcription cofactor vestigial like protein 4; Vestigial like protein 4; VGL 4; VGLL4; VGLL4_HUMAN

GENE INFORMATION

Entrez Gene ID [9686](#)

UniProt ID [Q14135](#)
