



Anti-PRSS41 (aa 105-133) polyclonal antibody (DPABH-02692)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PRSS41 (protease, serine, 41) is a protein-coding gene. Diseases associated with PRSS41 include focal epilepsy , and intellectual disability. GO annotations related to this gene include serine-type endopeptidase activity.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to internal sequence amino acids 105-133 of Human PRSS41 (NP_001128558.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name	PRSS41 protease, serine, 42 [Homo sapiens]
Official Symbol	PRSS41
Synonyms	PRSS41; protease, serine, 41; TESSP1; putative serine protease 41; testis-specific serine protease 1;
Entrez Gene ID	360226
Protein Refseq	NP_001128558.1
UniProt ID	Q7RTY9
Function	serine-type endopeptidase activity;