



Anti-SYTL5 (N-terminal) polyclonal antibody (DPABH-02795)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May act as Rab effector protein and play a role in vesicle trafficking. Binds phospholipids.
Specificity	DPABH-02795 is not predicted to cross-react with other SYTL protein family members.
Immunogen	Synthetic peptide corresponding to a 17 amino acid sequence near the N-terminus of Human SYTL5 (NP_620135).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC
Format	Liquid
Size	100 µg
Buffer	Constituent: 99.98% PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C, stable for one year. Should not be exposed to prolonged high temperatures. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SYTL5 synaptotagmin-like 6 [Homo sapiens]
Official Symbol	SYTL5
Synonyms	SYTL5; synaptotagmin-like 5; slp5; synaptotagmin-like protein 5; exophilin 9;
Entrez Gene ID	94122
Protein Refseq	NP_001156806.1
UniProt ID	Q8TDW5
Pathway	Deregulation of Rab and Rab Effector Genes in Bladder Cancer.
Function	Rab GTPase binding; metal ion binding; protein binding;