



Anti-TGFB3 (aa 250-350) polyclonal antibody (DPABH-02996)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Involved in embryogenesis and cell differentiation.
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 250 - 350 of Human TGF beta 3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituent: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	TGFB3 transforming growth factor, beta 4 [Homo sapiens]
Official Symbol	TGFB3
Synonyms	TGFB3; transforming growth factor, beta 3; ARVD; RNHF; ARVD1; TGF-beta3; transforming growth factor beta-3; TGF-beta-3; prepro-transforming growth factor beta-3;
Entrez Gene ID	7043
Protein Refseq	NP_003230.1
UniProt ID	A5YM40
Pathway	ALK1 signaling events; Amoebiasis; Cell cycle; Chagas disease (American trypanosomiasis)
Function	growth factor activity; identical protein binding; protein binding; contributes_to protein binding