



Anti-SH3GL1 (aa 261-272) polyclonal antibody (DPABH-03213)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Implicated in endocytosis. May recruit other proteins to membranes with high curvature.
Immunogen	Peptide with sequence C-EPFDLGEPEQSN, corresponding to internal sequence amino acids 261-272 of Human SH3 containing Grb 2 like 1 protein (NP_003016.1)
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	200 µl
Buffer	pH: 7.30; Constituents: 99% Tris buffered saline, 0.5% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [SH3GL1 SH3-domain GRB2-like 2 \[Homo sapiens \]](#)

Official Symbol	SH3GL1
Synonyms	SH3GL1; SH3-domain GRB2-like 1; EEN; CNSA1; SH3P8; SH3D2B; endophilin-A2; endophilin-2; SH3 domain protein 2B; EEN fusion partner of MLL; extra 11-19 leukemia fusion; SH3-containing Grb-2-like 1 protein; SH3 domain-containing GRB2-like protein 1; extra eleven-nineteen leukemia fusion gene protein;
Entrez Gene ID	6455
Protein Refseq	NP_001186872.1
UniProt ID	Q99961
Pathway	Endocytosis; Signaling events mediated by focal adhesion kinase;
Function	GTPase binding; identical protein binding; lipid binding; phosphatase binding