



Anti-RTTN (aa 50-150) polyclonal antibody (DPABH-03586)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Involved in the genetic cascade that governs left-right specification. Required for correct asymmetric expression of NODAL, LEFTY and PITX2.
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 50 - 150 of Human Rotatin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituent: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	RTTN rotatin [Homo sapiens]
Official Symbol	RTTN
Synonyms	RTTN; rotatin;
Entrez Gene ID	25914
Protein Refseq	NP_775901.3
UniProt ID	Q86VV8