



Anti-TGFB2 (aa 250-350) polyclonal antibody (DPABH-03686)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TGF-beta 2 has suppressive effects on interleukin-2 dependent T-cell growth.
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 250 - 350 of Human TGF beta 2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituent: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	TGFB2 transforming growth factor, beta 3 [Homo sapiens]
Official Symbol	TGFB2
Synonyms	TGFB2; transforming growth factor, beta 2; LDS4; TGF-beta2; transforming growth factor beta-2; G-TSF; cetermin; polyergin; BSC-1 cell growth inhibitor; prepro-transforming growth factor beta-2; glioblastoma-derived T-cell suppressor factor;
Entrez Gene ID	7042
Protein Refseq	NP_001129071.1
UniProt ID	P61812
Pathway	ATF-2 transcription factor network; Amoebiasis; Cell cycle; Chagas disease (American trypanosomiasis).
Function	beta-amyloid binding; cytokine activity; growth factor activity; protein binding