



Anti-TMEM70 polyclonal antibody (DPABH-03780)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Involved in biogenesis of mitochondrial ATP synthase.
Immunogen	Synthetic peptide of Human TMEM70, conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Format	Liquid
Size	50 µg
Buffer	Constituent: 99% PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name [TMEM70 transmembrane protein 71 \[Homo sapiens \]](#)

Official Symbol	TMEM70
Synonyms	TMEM70; transmembrane protein 70; MC5DN2; transmembrane protein 70, mitochondrial;
Entrez Gene ID	54968
Protein Refseq	NP_001035703.1
UniProt ID	Q9BUB7
Function	molecular_function;