



Anti-LARGE (aa 421-433) polyclonal antibody (DPABH-03888)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Glycosyltransferase which participates in glycosylation of alpha-dystroglycan. May carry out the synthesis of glycoprotein and glycosphingolipid sugar chains. May be involved in the addition of a repeated disaccharide unit.
Immunogen	Synthetic peptide: SEADVNSENLQKQ by a Cysteine residue linker, corresponding to internal sequence amino acids 421-433 of Human LARGE (NP_004728.1)
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	200 µl
Buffer	Constituents: 0.5% BSA, Tris buffered saline, pH 7.3
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	LARGE like-glycosyltransferase [Homo sapiens]
Official Symbol	LARGE
Synonyms	LARGE; like-glycosyltransferase; MDC1D; MDDGA6; MDDGB6; glycosyltransferase-like protein LARGE1; like-acetylglucosaminyltransferase; acetylglucosaminyltransferase-like 1A; acetylglucosaminyltransferase-like protein;
Entrez Gene ID	9215
Protein Refseq	NP_004728.1
UniProt ID	O95461
Function	acetylglucosaminyltransferase activity; transferase activity, transferring glycosyl groups;