



Anti-RGP1 (N-terminal) polyclonal antibody (DPABH-04029)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RGP1 (RGP1 retrograde golgi transport homolog (S. cerevisiae)) is a protein-coding gene. Diseases associated with RGP1 include periodontal disease , and periodontitis.
Immunogen	Synthetic peptide (17 amino acid peptide) near the N terminus of the Human RGP1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	Constituent: 99% PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name [RGP1 RGP2 retrograde golgi transport homolog \(S. cerevisiae\) \[Homo sapiens \]](#)

Official Symbol	RGP1
Synonyms	RGP1; RGP1 retrograde golgi transport homolog (S. cerevisiae); KIAA0258; retrograde Golgi transport protein RGP1 homolog;
Entrez Gene ID	9827
Protein Refseq	NP_001073965.2
UniProt ID	A8K0K1