



Anti-SEPT3 (aa 15-26) polyclonal antibody (DPABH-04094)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Filament-forming cytoskeletal GTPase (By similarity). May play a role in cytokinesis.
Immunogen	Synthetic peptide: C-SELVPEPRPKPA (with a Cysteine residue linker), corresponding to N terminal amino acids 15-26 of Human Septin 3
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	50 µg
Buffer	pH: 7.30; Constituents: 99% Tris buffered saline, 0.5% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name [SEPT3 septin 4 \[Homo sapiens \]](#)

Official Symbol	SEPT3
Synonyms	SEPT3; septin 3; SEP3; bK250D10.3; neuronal-specific septin-3; neuronal-specific septin 3;
Entrez Gene ID	55964
Protein Refseq	NP_061979.3
UniProt ID	Q9UH03
Pathway	Bacterial invasion of epithelial cells;
Function	GTP binding; molecular_function;