



Anti-OTULIN (aa 1-100) polyclonal antibody (DPABH-04254)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to Human FAM105B aa 1-100 conjugated to Keyhole Limpet Haemocyanin (KLH).Database link: Q96BN8
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituent: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [OTULIN OTU deubiquitinase with linear linkage specificity \[Homo sapiens \]](#)

Official Symbol	OTULIN
Synonyms	OTULIN; OTU deubiquitinase with linear linkage specificity; GUM; FAM105B; ubiquitin thioesterase otulin; gumby; protein FAM105B; ubiquitin thioesterase Gumby; deubiquitinating enzyme otulin; family with sequence similarity 105, member B; OTU domain-containing deubiquitinase with linear linkage specificity;
Entrez Gene ID	90268
Protein Refseq	NP_612357.4
UniProt ID	Q96BN8
Function	cysteine-type peptidase activity; protein binding; ubiquitin-specific protease activity;