



Anti-STXBP6 (aa 2-14) polyclonal antibody (DPABH-05010)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Forms non-fusogenic complexes with SNAP25 and STX1A and may thereby modulate the formation of functional SNARE complexes and exocytosis.
Immunogen	Synthetic peptide: SAKSAISKEIFAP-C, corresponding to N terminal amino acids 2-14 of Human Amisyn.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, IHC-P
Format	Liquid
Size	50 µg
Buffer	pH: 7.30; Constituents: 0.5% BSA, 99% Tris buffered saline
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	STXBP6 syntaxin binding protein 7 (amisyn) [Homo sapiens]
Official Symbol	STXBP6
Synonyms	STXBP6; syntaxin binding protein 6 (amisyn); amisyn; HSPC156; syntaxin-binding protein 6;
Entrez Gene ID	29091
Protein Refseq	NP_054897.4
UniProt ID	Q8NFX7