



Anti-FBXO9 (aa 431-447) polyclonal antibody (DPABH-05054)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.
Immunogen	Synthetic peptide: LFFARVRSYTAFSERPL conjugated to KLH, corresponding to C terminal amino acids 431-447 of Human FBXO9
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P, ICC/IF
Format	Liquid
Size	50 µg
Buffer	pH: 7.20; Constituents: 0.27% Potassium phosphate, 0.88% Sodium chloride
Preservative	0.01% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	FBXO9 F-box protein 10 [Homo sapiens]
Official Symbol	FBXO9
Synonyms	FBXO9; F-box protein 9; FBX9; VCIA1; NY-REN-57; dJ341E18.2; F-box only protein 9; F-box protein Fbx9; cross-immune reaction antigen 1; renal carcinoma antigen NY-REN-57;
Entrez Gene ID	26268
Protein Refseq	NP_036479.1
UniProt ID	Q9UK97
Function	protein binding; ubiquitin-protein transferase activity;