



# Rabbit Anti-Human XPO4 Polyclonal Antibody (DPABH-05184)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                                                                                                    |
|--------------------|--------------------------------------------------------------------------------------------------------------------|
| Immunogen          | Recombinant fusion protein containing a sequence corresponding to amino acids 400-700 of human XPO4 (NP_071904.4). |
| Isotype            | IgG                                                                                                                |
| Source/Host        | Rabbit                                                                                                             |
| Species Reactivity | Human                                                                                                              |
| Purification       | Affinity purification                                                                                              |
| Conjugate          | Unconjugated                                                                                                       |
| Applications       | WB                                                                                                                 |
| Format             | Liquid                                                                                                             |
| Size               | 50 µl, 100 µl                                                                                                      |
| Buffer             | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                                                                  |
| Preservative       | 0.02% Sodium Azide                                                                                                 |
| Storage            | Store at -20°C. Avoid freeze / thaw cycles.                                                                        |

## BACKGROUND

|          |                                     |
|----------|-------------------------------------|
| Keywords | XPO4; exportin 4; exp4; exportin-4; |
|----------|-------------------------------------|

# GENE INFORMATION

Entrez Gene ID [64328](#)

---

UniProt ID [Q9C0E2](#)

---