



Anti-SLITRK2 (C-terminal) polyclonal antibody (DPABH-05196)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Suppresses neurite outgrowth.
Immunogen	A 17 amino acid synthetic peptide from near the carboxy terminus of Human SLITRK2 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	50 µg
Buffer	Constituent: 99% PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name [SLITRK2 SLIT and NTRK-like family, member 3 \[Homo sapiens \]](#)

Official Symbol	SLITRK2
Synonyms	SLITRK2; SLIT and NTRK-like family, member 2; CXorf2; SLITL1; SLIT and NTRK-like protein 2; slit-like 1; slit and trk like gene 2;
Entrez Gene ID	84631
Protein Refseq	NP_001137475.1
UniProt ID	Q9H156