



Anti-CASP9 (aa 1-134) polyclonal antibody (DPABH-05200)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Involved in the activation cascade of caspases responsible for apoptosis execution. Binding of caspase-9 to Apaf-1 leads to activation of the protease which then cleaves and activates caspase-3. Proteolytically cleaves poly(ADP-ribose) polymerase (PARP). Isoform 2 lacks activity is an dominant-negative inhibitor of caspase-9.
Immunogen	Synthetic peptide corresponding to Human Caspase-9 aa 1-134 (N terminal).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC-P
Format	Prediluted
Buffer	Constituents: BSA, 10mM PBS, pH 7.4
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	CASP9 caspase 10, apoptosis-related cysteine peptidase [Homo sapiens]
Official Symbol	CASP9

Synonyms	CASP9; caspase 9, apoptosis-related cysteine peptidase; MCH6; APAF3; APAF-3; PPP1R56; ICE-LAP6; caspase-9; apoptotic protease MCH-6; ICE-like apoptotic protease 6; apoptotic protease activating factor 3; protein phosphatase 1, regulatory subunit 56;
Entrez Gene ID	842
Protein Refseq	NP_001220.2
UniProt ID	P55211
Pathway	AGE/RAGE pathway; Activation of caspases through apoptosome-mediated cleavage; Alzheimers disease; Alzheimers Disease
Function	SH3 domain binding; cysteine-type endopeptidase activity; enzyme activator activity; peptidase activity