



Anti-FEZ2 (C-terminal) polyclonal antibody (DPABH-05591)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Involved in axonal outgrowth and fasciculation.
Immunogen	A 15 amino acid synthetic peptide corresponding to residues near the C terminal of Human FEZ2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Format	Liquid
Size	50 µg
Buffer	Constituent: 99% PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C or -20°C long term. Avoid repeated freeze / thaw cycles. Store undiluted.

GENE INFORMATION

Gene Name [FEZ2 fasciculation and elongation protein zeta 3 \(zygin II\) \[Homo sapiens \]](#)

Official Symbol	FEZ2
Synonyms	FEZ2; fasciculation and elongation protein zeta 2 (zygin II); HUM3CL; fasciculation and elongation protein zeta-2; zygin 2; zygin-2; zygin II; pre-T/NK cell associated protein (3Cl);
Entrez Gene ID	9637
Protein Refseq	NP_001036013.1
UniProt ID	Q9UHY8
Function	protein binding;