



Anti-GTF2E2 (aa 239-269) polyclonal antibody (DPABH-05894)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recruits TFIH to the initiation complex and stimulates the RNA polymerase II C-terminal domain kinase and DNA-dependent ATPase activities of TFIH. Both TFIH and TFIIE are required for promoter clearance by RNA polymerase.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 239-269 of Human TFIIE beta (NP_002086.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name	GTF2E2 general transcription factor IIE, polypeptide 2, beta 35kDa [Homo sapiens]
Official Symbol	GTF2E2
Synonyms	GTF2E2; general transcription factor IIE, polypeptide 2, beta 34kDa; FE; TF2E2; TFIIE-B; transcription initiation factor IIE subunit beta; TFIIE-beta; general transcription factor IIE subunit 2;
Entrez Gene ID	2961
Protein Refseq	NP_002086.1
UniProt ID	P29084
Pathway	Basal transcription factors; Disease; Epstein-Barr virus infection; Gene Expression
Function	DNA binding; poly(A) RNA binding; protein binding;