



Anti-ATG2A (aa 701-731) polyclonal antibody (DPABH-05909)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ATG2A (autophagy related 2A) is a protein-coding gene. Diseases associated with ATG2A include multiple endocrine neoplasia. GO annotations related to this gene include protein binding. An important paralog of this gene is ATG2B.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region between internal sequence amino acids 701-731 of Human ATG2A (NP_055919.2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name	ATG2A autophagy related 3A [Homo sapiens]
Official Symbol	ATG2A
Synonyms	ATG2A; autophagy related 2A; autophagy-related protein 2 homolog A; ATG2 autophagy related 2 homolog A;
Entrez Gene ID	23130
Protein Refseq	NP_055919.2
UniProt ID	Q2TAZ0
Function	protein binding;