



Anti-ANKS6 (aa 504-534) polyclonal antibody (DPABH-05925)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Required for renal function.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 504-534 of Human ANKS6 (NP_775822.3).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name [ANKS6 ankyrin repeat and sterile alpha motif domain containing 7 \[Homo sapiens \]](#)

Official Symbol	ANKS6
Synonyms	ANKS6; ankyrin repeat and sterile alpha motif domain containing 6; PKDR1; SAMD6; NPHP16; ANKRD14; ankyrin repeat and SAM domain-containing protein 6; samCystin; ankyrin repeat domain 14; SAM domain-containing protein 6; sterile alpha motif domain containing 6; ankyrin repeat domain-containing protein 14; sterile alpha motif domain-containing protein 6;
Entrez Gene ID	203286
Protein Refseq	NP_775822.3
UniProt ID	B3KXP1
Function	protein homodimerization activity;