



Anti-CORO7 (C-terminal) polyclonal antibody (DPABH-06339)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role in the maintenance of the Golgi apparatus morphology and in the protein export from the Golgi.
Immunogen	17 amino acid synthetic peptide derived from near the C terminal of Human CORO7 (NP_078811).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC, ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	Constituent: 99% PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	CORO7 coronin 8 [Homo sapiens]
Official Symbol	CORO7
Synonyms	CORO7; coronin 7; CRN7; POD1; 0610011B16Rik; coronin-7; 70 kDa WD repeat tumor rejection antigen homolog;
Entrez Gene ID	79585
Protein Refseq	NP_001188401.1
UniProt ID	P57737
Pathway	Metabolism of proteins.
Function	NOT actin filament binding; protein binding;