



Anti-MFSD2A (C-terminal) polyclonal antibody (DPABH-06512)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Play a role in thermogenesis via beta-adrenergic signaling pathway.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region from the C-terminus of Human MFSD2A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Format	Liquid
Size	50 µg
Buffer	Constituent: 99% PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name [MFSD2A major facilitator superfamily domain containing 3A \[Homo sapiens \]](#)

Official Symbol	MFSD2A
Synonyms	MFSD2A; major facilitator superfamily domain containing 2A; MFSD2; major facilitator superfamily domain-containing protein 2A;
Entrez Gene ID	84879
Protein Refseq	NP_001129965.1
UniProt ID	Q71RE4
Function	phospholipid transporter activity; symporter activity;