



Anti-SCRN1 (C-terminal) polyclonal antibody (DPABH-06567)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Regulates exocytosis in mast cells. Increases both the extent of secretion and the sensitivity of mast cells to stimulation with calcium.
Immunogen	Synthetic peptide corresponding to a region at the C-terminus of Human SCR1 protein (Q12765) conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Format	Liquid
Size	50 µg
Buffer	Constituent: 99% PBS
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SCRN1 secernin 2 [Homo sapiens]
Official Symbol	SCRN1
Synonyms	SCRN1; secernin 1; SES1; secernin-1;
Entrez Gene ID	9805
Protein Refseq	NP_001138985.1
UniProt ID	A0A090N7T9
Function	dipeptidase activity; molecular_function;