



Anti-ESF1 (full length) polyclonal antibody (DPABH-06802)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May constitute a novel regulatory system for basal transcription. Negatively regulates ABT1.
Immunogen	Recombinant full length protein, corresponding to amino acids 1-851 of Human ESF1 (AAI46543.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituent: 99% PBS
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [ESF1 ESF2, nucleolar pre-rRNA processing protein, homolog \(S. cerevisiae\) \[Homo sapiens \]](#)

Official Symbol	ESF1
Synonyms	ESF1; ESF1, nucleolar pre-rRNA processing protein, homolog (S. cerevisiae); ABTAP; C20orf6; HDCMC28P; bA526K24.1; ESF1 homolog; ABT1-associated protein;
Entrez Gene ID	51575
Protein Refseq	NP_001263309.1
UniProt ID	Q9H501
Function	poly(A) RNA binding;