



Anti-G0S2 (full length) polyclonal antibody (DPABH-06934)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Potential oncogene and regulator of latent HIV.
Immunogen	Recombinant full length protein, corresponding to amino acids 1-103 of Human G0/G1switch 2 (AAH09694.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	G0S2 G0/G1 switch 3 [Homo sapiens]
Official Symbol	G0S2
Synonyms	G0S2; G0/G1 switch 2; G0/G1 switch protein 2; G0/G1switch 2; G0/G1 switch regulatory protein 2;
Entrez Gene ID	50486
Protein Refseq	NP_056529.1
UniProt ID	P27469
Pathway	Diurnally regulated genes with circadian orthologs; Metabolism; PPARA activates gene expression;
Function	molecular_function; protein binding;