



Anti-KIR2DL5A (full length) polyclonal antibody (DPABH-06979)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Receptor on natural killer (NK) cells for HLA-C alleles. Inhibits the activity of NK cells thus preventing cell lysis.
Immunogen	Recombinant full length protein, corresponding to amino acids 1-375 of Human KIR2DL5A (AAI60063.1).
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	pH: 7.20; Constituent: 99% PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	KIR2DL5A killer cell immunoglobulin-like receptor, two domains, long cytoplasmic tail, 6A [Homo sapiens]
Official Symbol	KIR2DL5A
Synonyms	KIR2DL5A; killer cell immunoglobulin-like receptor, two domains, long cytoplasmic tail, 5A; CD158F; KIR2DL5; KIR2DL5.1; KIR2DL5.3; killer cell immunoglobulin-like receptor 2DL5A; killer-cell Ig-like receptor; killer cell immunoglobulin-like receptor KIR2DL5A; killer cell immunoglobulin-like receptor, two domains, long cytoplasmic tail, 5;
Entrez Gene ID	57292
Protein Refseq	NP_065396.1
UniProt ID	Q8N109
Pathway	Antigen processing and presentation; Graft-versus-host disease; Natural killer cell mediated cytotoxicity;