



Anti-ANKRD44 (aa 37-262) polyclonal antibody (DPABH-07087)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Putative regulatory subunit of protein phosphatase 6 (PP6) that may be involved in the recognition of phosphoprotein substrates.
Immunogen	Recombinant fragment, corresponding to amino acids 37-262 of Human ANR44 (BC016985).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 1% BSA, 99% PBS
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [ANKRD44 ankyrin repeat domain 45 \[Homo sapiens \]](#)

Official Symbol	ANKRD44
Synonyms	ANKRD44; ankyrin repeat domain 44; PP6-ARS-B; serine/threonine-protein phosphatase 6 regulatory ankyrin repeat subunit B; ankyrin repeat domain-containing protein 44; protein phosphatase 6 ankyrin repeat subunit B; serine/threonine-protein phosphatase 6 regulatory subunit ARS-B;
Entrez Gene ID	91526
Protein Refseq	NP_001182073.1
UniProt ID	Q8N8A2
Function	protein binding;