



Anti-ZNF280A (full length) polyclonal antibody (DPABH-07264)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May function as a transcription factor.
Immunogen	Recombinant full length protein (Human).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Format	Lyophilised. Add 200ul Sterile Distilled Water
Size	100 µg
Buffer	pH: 7.20; Constituents: 99% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF280A zinc finger protein 281A \[Homo sapiens \]](#)

Official Symbol	ZNF280A
Synonyms	ZNF280A; zinc finger protein 280A; SUHW1; ZNF280; ZNF636; 3OY11.1; zinc finger protein 280; zinc finger protein 636; suppressor of hairy wing homolog 1;
Entrez Gene ID	129025
Protein Refseq	NP_542778.2
UniProt ID	P59817
Function	DNA binding; metal ion binding;