



Rabbit anti-Human ANKRD40 polyclonal antibody (DPABH-07353)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ANKRD40 (ankyrin repeat domain 40) is a protein-coding gene. GO annotations related to this gene include molecular_function.
Immunogen	Recombinant fragment, corresponding to amino acids 126-326 of Human ANKRD40 (BC012978).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 98% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ANKRD40 ankyrin repeat domain 41 [Homo sapiens]
Official Symbol	ANKRD40
Synonyms	ANKRD40; ankyrin repeat domain 40; ankyrin repeat domain-containing protein 40;
Entrez Gene ID	91369
Protein Refseq	NP_443087.1
UniProt ID	A8IK34
Function	molecular_function;