



Rabbit Anti-Human ANKDD1A Polyclonal Antibody (DPABH-07364)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: LTALHSAAGGSHPCVQLLLRAGSTVNALTQKNLSCLHYAALSGSEDVSRVLIHAGGCAN VVDHQGASPLHLAVRHNFPALVRLINSDSDVNAV
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purified
Conjugate	Unconjugated
Applications	IHC, ICC-IF
Format	Liquid
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2).
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

BACKGROUND

Keywords

ANKDD1A; ankyrin repeat and death domain containing 1A; ankyrin repeat and death domain-containing protein 1A;

GENE INFORMATION

Entrez Gene ID[348094](#)**UniProt ID**[Q495B1](#)