



Anti-JOSD1 (aa 44-188) polyclonal antibody (DPABH-07419)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Has low protease activity towards poly-ubiquitin chains (in vitro). May act as a deubiquitinating enzyme.
Immunogen	Recombinant fragment, corresponding to amino acids 44-188 of Human Josephin-1 (BC015026).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 98% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	JOSD1 Josephin domain containing 2 [Homo sapiens]
Official Symbol	JOSD1
Synonyms	JOSD1; Josephin domain containing 1; dJ508I15.2; Josephin-1; josephin domain-containing 1; josephin domain-containing protein 1;
Entrez Gene ID	9929
Protein Refseq	NP_055691.1
UniProt ID	A0A024R1P5
Function	omega peptidase activity;