



Anti-ZNF404 (aa 63-314) polyclonal antibody (DPABH-07551)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	Recombinant fragment, corresponding to amino acids 63-314 of Human ZNF404 (BC101333).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 98% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF404 zinc finger protein 405 \[Homo sapiens \]](#)

Official Symbol	ZNF404
Synonyms	ZNF404; zinc finger protein 404;
Entrez Gene ID	342908
Protein Refseq	NP_001028891.2
UniProt ID	Q494X3
Function	DNA binding; metal ion binding;