



# Anti-ZNF512 (aa 16-268) polyclonal antibody (DPABH-07562)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Antigen Description | May be involved in transcriptional regulation.                                                 |
| Immunogen           | Recombinant fragment, corresponding to amino acids 16-268 of Human ZNF512 (BC043221).          |
| Isotype             | IgG                                                                                            |
| Source/Host         | Rabbit                                                                                         |
| Species Reactivity  | Human                                                                                          |
| Purification        | Immunogen affinity purified                                                                    |
| Conjugate           | Unconjugated                                                                                   |
| Applications        | WB                                                                                             |
| Format              | Liquid                                                                                         |
| Size                | 100 µg                                                                                         |
| Buffer              | pH: 7.20; Constituents: 1% BSA, 98% PBS                                                        |
| Preservative        | 0.02% Sodium Azide                                                                             |
| Storage             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |

## GENE INFORMATION

Gene Name [ZNF512 zinc finger protein 513 \[ Homo sapiens \]](#)

|                        |                                  |
|------------------------|----------------------------------|
| <b>Official Symbol</b> | ZNF512                           |
| <b>Synonyms</b>        | ZNF512; zinc finger protein 512; |
| <b>Entrez Gene ID</b>  | <a href="#">84450</a>            |
| <b>Protein Refseq</b>  | <a href="#">NP_001258215.1</a>   |
| <b>UniProt ID</b>      | <a href="#">Q96ME7</a>           |
| <b>Function</b>        | DNA binding; metal ion binding;  |