



Anti-ZNF322 (aa 139-402) polyclonal antibody (DPABH-07573)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Involved in transcriptional regulation as an activator.
Immunogen	Recombinant protein Human ZNF322A ,corresponding to amino acids 139-402 of Human ZNF322A
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 98% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF322 zinc finger protein 323 \[Homo sapiens \]](#)

Official Symbol	ZNF322
Synonyms	ZNF322; zinc finger protein 322; HCG12; ZNF388; ZNF489; ZNF322A; HLA complex group 12; zinc finger protein 388; zinc finger protein 489; zinc finger protein 322A;
Entrez Gene ID	79692
Protein Refseq	NP_001229726.1
UniProt ID	A0A024QZZ4
Function	DNA binding; metal ion binding; sequence-specific DNA binding transcription factor activity;