



Anti-SCARB2 (aa 25-74) polyclonal antibody (DPABH-08218)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Acts as a lysosomal receptor for glucosylceramidase (GBA) targeting.
Immunogen	Synthetic peptide corresponding a region within N terminal amino acids 25-74 (VARVFQKAVD QSIEKKIVLR NGTEAFDSWE KPPLPVYTQF YFFNVTNPEE) of Human LIMP2 (NP_005497).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [SCARB2 scavenger receptor class B, member 3 \[Homo sapiens \]](#)

Official Symbol	SCARB2
Synonyms	SCARB2; scavenger receptor class B, member 2; AMRF; EPM4; LGP85; CD36L2; HLGP85; LIMP-2; LIMP2; SR-BII; lysosome membrane protein 2; LIMP II; CD36 antigen-like 2; lysosome membrane protein II; 85 kDa lysosomal membrane sialoglycoprotein; 85 kDa lysosomal sialoglycoprotein scavenger receptor class B, member 2; CD36 antigen (collagen type I receptor, thrombospondin receptor)-like 2 (lysosomal integral membrane protein II);
Entrez Gene ID	950
Protein Refseq	NP_001191184.1
UniProt ID	Q14108
Pathway	Lysosome;
Function	enzyme binding; protein binding; receptor activity;