



Rabbit Anti-Human POU2AF1 Polyclonal Antibody (DPABH-08449)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal to Bob 1.
Specificity	The antibody detects endogenous Bob 1.
Target	Bob 1
Immunogen	Synthetic peptide from human Bob 1 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	This antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Conjugate	Unconjugated
Applications	IHC-P, ELISA
Format	Liquid
Concentration	Lot specific
Size	100 µl, 200 µl
Buffer	PBS, pH 7.4, containing 50% Glycerol.
Preservative	0.02% Sodium Azide

Storage

Store at -20°C, and avoid repeat freeze-thaw cycles.

BACKGROUND

Introduction

POU2AF1 (POU class 2 associating factor 1) is a protein-coding gene. Diseases associated with POU2AF1 include primary biliary cirrhosis, and primary effusion lymphoma. GO annotations related to this gene include transcription coactivator activity and transcription cofactor activity.

Keywords

POU2AF1;POU class 2 associating factor 1;BOB1;OBF1;OCAB;OBF-1;POU domain class 2-associating factor 1;BOB-1;OCA-B;OCT-binding factor 1;B-cell-specific coactivator OBF-1;POU domain, class 2, associating factor 1;

GENE INFORMATION

Entrez Gene ID

[5450](#)

UniProt ID

[Q16633](#)
