



Rabbit Anti-Human PAPD5 Polyclonal Antibody (DPABH-08492)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: QSSSSDVSDATPCKTPKQLLCRPSTGNRVGSQDVSLESSQAVGKMQSTQTTNTSNSTNK SQHGSARLFRSSSKGFQGTQTSH
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2).
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

BACKGROUND

Introduction

PAPD5 (PAP associated domain containing 5) is a protein-coding gene. GO annotations related to this gene include DNA-directed DNA polymerase activity and molecular_function. An important paralog of this gene is PAPD7.

Keywords

PAPD5; PAP associated domain containing 5; TRF4-2; non-canonical poly(A) RNA polymerase PAPD5; TUTase 3; terminal uridylyltransferase 3; PAP-associated domain-containing protein 5; topoisomerase-related function protein 4-2;

GENE INFORMATION

Entrez Gene ID

[64282](#)

UniProt ID

[Q8NDF8](#)