



Anti-C11ORF95 (aa 563-592) polyclonal antibody (DPABH-08551)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide within Human C11orf95 aa 563-592 (C terminal) conjugated to Keyhole Limpet Haemocyanin (KLH). The exact sequence is proprietary.Database link: C9JLR9
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB, ICC/IF
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name [C11ORF95 chromosome 11 open reading frame 96 \[Homo sapiens \]](#)

Official Symbol	C11ORF95
Synonyms	C11ORF95; chromosome 11 open reading frame 95; uncharacterized protein C11orf95;
Entrez Gene ID	65998
Protein Refseq	NP_001138408.1
UniProt ID	C9JLR9