



Anti-OR4F15 (aa 76-105) polyclonal antibody (DPABH-08554)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Odorant receptor.
Immunogen	Synthetic peptide within Human O4F15 aa 76-105 conjugated to Keyhole Limpet Haemocyanin (KLH). The exact sequence is proprietary.Database link: Q8NGB8
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, Flow Cyt
Format	Liquid
Size	100 µl
Buffer	Constituent: 99.91% PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	OR4F15 olfactory receptor, family 4, subfamily F, member 16 [Homo sapiens]
Official Symbol	OR4F15
Synonyms	OR4F15; olfactory receptor, family 4, subfamily F, member 15; olfactory receptor 4F15; olfactory receptor OR15-14;
Entrez Gene ID	390649
Protein Refseq	NP_001001674.1
UniProt ID	Q8NGB8
Pathway	GPCR downstream signaling; Olfactory transduction; Signal Transduction;
Function	G-protein coupled receptor activity; olfactory receptor activity;