



Anti-KHDC3L (aa 27-56) polyclonal antibody (DPABH-08563)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KHDC3L (KH domain containing 3-like, subcortical maternal complex member) is a protein-coding gene. Diseases associated with KHDC3L include hydatidiform mole, recurrent. An important paralog of this gene is OOEP.
Immunogen	Synthetic peptide within Human Ecat1 aa 27-56 (N terminal) conjugated to Keyhole Limpet Haemocyanin (KLH). The exact sequence is proprietary. Database link: Q587J8
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, Flow Cyt, IHC-P
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	KHDC3L KH domain containing 4-like, subcortical maternal complex member [Homo sapiens]
Official Symbol	KHDC3L
Synonyms	KHDC3L; KH domain containing 3-like, subcortical maternal complex member; ECAT1; HYDM2; C6orf221; KHDC3-like protein; ES cell-associated transcript 1 protein;
Entrez Gene ID	154288
Protein Refseq	NP_001017361.1
UniProt ID	Q587J8
Function	RNA binding;