



Anti-ZKSCAN4 (aa 183-476) polyclonal antibody (DPABH-08602)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in the transcriptional activation of MDM2 and EP300 genes.
Immunogen	Recombinant fragment, corresponding to amino acids 183-476 of Human ZNF307 (BC014031).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 98% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZKSCAN4 zinc finger with KRAB and SCAN domains 5 \[Homo sapiens \]](#)

Official Symbol	ZKSCAN4
Synonyms	ZKSCAN4; zinc finger with KRAB and SCAN domains 4; ZNF307; ZNF427; ZSCAN36; P1P373C6; p373c6.1; zinc finger protein with KRAB and SCAN domains 4; zinc finger protein 307; zinc finger protein 427; p373c6.1 (novel C2H2 type zinc finger protein);
Entrez Gene ID	387032
Protein Refseq	NP_061983.2
UniProt ID	Q969J2
Pathway	Gene Expression;
Function	DNA binding; metal ion binding; sequence-specific DNA binding transcription factor activity;