



Anti-PLCL2 (aa 1100-1127) polyclonal antibody (DPABH-08911)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play an role in the regulation of Ins(1,4,5)P3 around the endoplasmic reticulum.
Immunogen	Synthetic peptide within Human PLCL2 aa 1100-1127 (C terminal) conjugated to Keyhole Limpet Haemocyanin (KLH). The exact sequence is proprietary.Database link: Q9UPRO
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P, Flow Cyt, ICC/IF
Format	Liquid
Size	100 µg
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	PLCL2 phospholipase C-like 2 [Homo sapiens]
Official Symbol	PLCL2
Synonyms	PLCL2; phospholipase C-like 2; phospholipase C, epsilon 2; PLCE2; inactive phospholipase C-like protein 2; KIAA1092; PLC-L2; PLC-L(2); PLC-epsilon-2; phospholipase C-L2; phospholipase C-epsilon-2; phospholipase C, epsilon 2; PLCE2; FLJ13484;
Entrez Gene ID	23228
Protein Refseq	NP_001137854
UniProt ID	Q9UPR0
Chromosome Location	3p25.3-p25.1
Function	calcium ion binding; phosphatidylinositol phospholipase C activity; signal transducer activity;