



Anti-ZNF20 (aa 249-522) polyclonal antibody (DPABH-08984)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	Recombinant fragment, corresponding to amino acids 249-522 of Human ZNF20 (BC036714).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 98% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF20 zinc finger protein 21 \[Homo sapiens \]](#)

Official Symbol	ZNF20
Synonyms	ZNF20; zinc finger protein 20; KOX13; zinc finger protein KOX13;
Entrez Gene ID	7568
Protein Refseq	NP_001190179.1
UniProt ID	P17024
Pathway	Gene Expression;
Function	DNA binding; metal ion binding;